



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

Table of contents

Generalities.....	2
1. Objectives.....	2
2. References	2
3. Definitions, Terms and Abbreviations.....	2
4. Signature table	3
Installation qualification (IQ)	4
5. Objective.....	4
6. Scope	4
7. Instrument origin and identification.....	4
8. Instrument delivery	5
9. Environment and siting	6
10. Instrument safety	6
11. Installation	7
12. IQ summary report.....	8
Operational Qualification (OQ)	9
13. Objective.....	9
14. Scope	9
15. Functional tests	10
16. Training.....	11
17. OQ summary report.....	12



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

Generalities

1. Objectives

The objectives of the IQ are to ensure that equipment documentation, delivery and installation conform to the manufacturer's literature, and to confirm the reference information of the equipment.

The objectives of the OQ are to determine that the equipment operates according to specifications, and to record all relevant information and data to demonstrate it functions as expected according to the manufacturer's specifications.

When IQ/OQ protocol is approved after execution, the instrument will be ready for use in non-critical applications and for further testing.

2. References

2.1 The following documents are referred to within the IQ/OQ protocol

No	Document title	Document ID	Comments
1	Miris HMA™ User Manual	DOC-510/DOC-511	
2	Miris Calibration Control Kit™ Instructions for Use	DOC-551/DOC-552	

3. Definitions, Terms and Abbreviations

3.1 Definitions

Abbreviation	Description
CPG	Cuvette Pressure Guard
HMA	Human Milk Analyzer
IFU	Instructions for Use
IQ/OQ	Installation Qualification / Operational Qualification



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

4. Signature table

All the trained personnel who have been working on the protocol must fill the information required in the list below and put down their signatures.

5.1 Signatures

Name	Title	Signature	Trained according to [1]	Initials



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

Installation qualification (IQ)

5. Objective

The objectives of the IQ are to ensure that equipment documentation, delivery and installation conform to the manufacturer's literature, and to confirm the reference information of the equipment.

6. Scope

IQ should be performed at the time of installation, modification, or relocation. The instrument shall be in power-off during the execution of the IQ protocol.

7. Instrument origin and identification

7.1 Instrument references

Instrument name	Miris HMA™ - Human Milk Analyzer
Instrument description	Mid-infrared transmission spectroscopy instrument for determining the macronutrient composition (fat, protein, carbohydrate, total solids, energy) of human milk
Manufacturer	Miris AB Danmarksgatan 26 753 23 UPPSALA SWEDEN www.MirisSolutions.com
Service contact	Tel: +46 18 14 69 07 e-mail: support@MirisSolutions.com
Instrument serial number product number software version	
Date received	
Condition when received	New <input type="checkbox"/> Serviced <input type="checkbox"/>
Location installed	
<i>Verified by:</i> <i>Date and signature:</i>	



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

8. Instrument delivery

Unpack and check the contents of the instrument carrying case. Please note that when an instrument is returned after a repair or service, all items listed here may not be included.

8.1 Are any of the standard items missing?

Item	Present	Missing
Miris HMA™ - Human Milk Analyzer	<input type="checkbox"/>	<input type="checkbox"/>
AC/DC Adapter	<input type="checkbox"/>	<input type="checkbox"/>
Power cable	<input type="checkbox"/>	<input type="checkbox"/>
Wireless mouse including batteries	<input type="checkbox"/>	<input type="checkbox"/>
USB flash drive	<input type="checkbox"/>	<input type="checkbox"/>
USB hub	<input type="checkbox"/>	<input type="checkbox"/>
Outlet Tubes, 6 pcs	<input type="checkbox"/>	<input type="checkbox"/>
Spare parts CPG	<input type="checkbox"/>	<input type="checkbox"/>
Miris HMA™ User Manual	<input type="checkbox"/>	<input type="checkbox"/>

Manufacturer informed of any missing items	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Comments and follow-up		
Verified by: Date and signature:		

8.2 Are there any damages to the instrument or accessories?

No <input type="checkbox"/>	Yes <input type="checkbox"/>	
	If yes, name the item and describe the damage	
	Corrective action	
	Manufacturer informed	Yes <input type="checkbox"/>
Verified by: Date and signature:		



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

9. Environment and siting

9.1 Since the instrument was last subjected to a change in environmental conditions, has the instrument been allowed to adjust to room temperature (20 - 30°C/68 - 86°F) for a minimum of 4 hours?

Yes No If no, ensure this condition is met

9.2 Verify that the operating environment of the instrument meets the below requirements

Operating parameter	Specified range	Condition met	Deviation
Ambient temperature	20 - 30° C (68 - 86° F)	<input type="checkbox"/>	
Humidity	20 - 80% not condensed	<input type="checkbox"/>	
Power supply adapter	Input voltage 100-240 V, ~ 50/60 Hz, 2.3A	<input type="checkbox"/>	
Power supply instrument	Output voltage 18 V DC, 100VA	<input type="checkbox"/>	
<i>Comments</i>			
<i>Verified by: Date and signature:</i>			

10. Instrument safety

10.1 Verify that the Manufacturer's safety recommendations are met

Condition	Verified	Comment
Allow enough space for the instrument	<input type="checkbox"/>	
Place on a level, stable surface	<input type="checkbox"/>	
Place in an area free from dust, dirt, explosives, corrosive fumes, and extremes of temperature and humidity	<input type="checkbox"/>	
Place to avoid draft and vibrations	<input type="checkbox"/>	
Do not place in direct sunlight	<input type="checkbox"/>	
Ensure correct power supply	<input type="checkbox"/>	
In case of unstable power supply, use an UPS (uninterruptible power supply)	<input type="checkbox"/>	N/A <input type="checkbox"/>
Leave the protective cap covering the RS232 connection when not in use	<input type="checkbox"/>	
Read the User Manual before use	<input type="checkbox"/>	

<i>Further comments and follow-up</i>	
<i>Verified by: Date and signature:</i>	



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

11. Installation

11.1 The installation is performed by: user specialised technician/engineer

11.2 Perform/verify the following installation procedure

Item	OK	Deviation
Unpack the instrument (4. Instrument delivery and documentation)	<input type="checkbox"/>	<input type="checkbox"/>
Control of the environment and instrument placement (5. Environment and siting, 6. Instrument safety)	<input type="checkbox"/>	<input type="checkbox"/>
Connect the instrument to power supply (see User Manual Chapter 1, Electrical requirements)	<input type="checkbox"/>	<input type="checkbox"/>
Attach tubes to the outlet (see User Manual Chapter 1, Operating the instrument)	<input type="checkbox"/>	<input type="checkbox"/>
Place a waste container by the outlet	<input type="checkbox"/>	<input type="checkbox"/>

<i>Comments including any deviations</i>	
<i>Verified by:</i>	
<i>Date and signature:</i>	



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

12. IQ summary report

Instrument: Miris HMA™ - Human Milk Analyzer

Manufacturer: Miris AB

Instrument serial number:

12.1 Assessment of complete Installation Qualification:

No deviations

Deviations (list below)

12.2 IQ deviations

Deviation No	Deviation	Impact on operation	Justification for acceptance

Completion of the preceding protocol indicates that this instrument has been adequately delivered and installed. The Installation Qualification has been passed and the instrument may be submitted for Operational Qualification.

Deviations approved by:		Date:	
IQ approved by:		Date:	



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

Operational Qualification (OQ)

13. Objective

Objectives of the OQ is to determine that the equipment operates according to specifications, and to record all relevant information and data to demonstrate it functions as expected according to the manufacturer's specifications.

14. Scope

OQ should be performed after installation, modification or relocation, and after the Installation Qualification has been approved.



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

15. Functional tests

15.1 Perform the following tests and checks sequentially

Action	Refer to	Acceptance criteria	Pass	Fail
Switch on the Miris HMA™ (ON/OFF switch)	-	Display illuminates and button indicator lights comes on	<input type="checkbox"/>	<input type="checkbox"/>
Fan (no action required)	-	Fan comes on at power-up	<input type="checkbox"/>	<input type="checkbox"/>
System load (no action required)	User Manual Chapter 1, The start-up procedure	1. Main menu appears 2. Message "System loaded" is displayed	<input type="checkbox"/>	<input type="checkbox"/>
Thermal regulation (no action required)	User Manual Chapter 1, The start-up procedure	1. Message "Warming..." is displayed 2. Temperature indicator at top right corner of display reaches 40°C 3. Message "Ready - press a button" is displayed	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Perform the instrument start-up procedure: 1. 30 minutes switched on 2. Inject cleaning solution 3. Inject distilled or deionized water	User Manual Chapter 1, The start-up procedure	N/A	Performed <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
Calibration load (no action required)	-	HMA Homogenized milk (1) is displayed at the top of the display	<input type="checkbox"/>	<input type="checkbox"/>
Zero-setting check (zero level)	User Manual Chapter 4, Check procedure with zero level adjustment	1. Change% within ±10% 2. Result% within ±0.05 3. Message "No adjustment necessary" is displayed	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Control of the calibration	IFU Miris Calibration Control Kit™	Measurement results for the Calibration 1 setting are within the acceptance criteria, for: Fat Crude protein True protein Carbohydrate	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
<i>Comments</i>				
<i>Verified by: Date and signature:</i>				

If any of the functional tests fails, contact Miris AB or your distributor.



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

16. Training

16.1 Training verification

Trained operators	Read the User Manual	Watched instructional videos (www.MirisSolutions.com)	Practical training	Authorising signature
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<i>Comments</i>				
<i>Verified by: Date and signature:</i>				



IQ OQ		
	Document title	
	Installation Qualification and Operational Qualification Miris HMA	

17. OQ summary report

Instrument: Miris HMA™ - Human Milk Analyzer

Manufacturer: Miris AB

Instrument serial number:

17.1 Assessment of complete Operational Qualification:

No deviations

Deviations (list below)

17.2 OQ deviations

Deviation No	Deviation	Impact on operation	Justification for acceptance

Completion of the preceding activities and checks indicates that this instrument is operating satisfactorily following delivery and installation. The instrument has passed the Operational Qualification procedure and may now be released for use in non-critical applications and further testing.

Deviations approved by:		Date:	
OQ approved by:		Date:	